Patella Tracking Orthosis (PTO) Knee Brace

PTO-KNIESCHIENE

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


E/U authorized representative
MDSS
Burckhardtstrasse 1
D-30163 Hanover
Germany

BREG Inc.
Vista, CA 92081 U.S.A.
Telephone (800) 321-0607
(760) 599-3000
Telefax (760) 598-6193

PN-1.08770 REV E

WARNING

CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PRACTITIONER.

CAUTION: FOR SINGLE PATIENT USE ONLY.

CAUTION: IF YOU EXPERIENCE INCREASED PAIN, SWELLING, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.

WARNING: CAREFULLY READ USE INSTRUCTIONS AND WARNINGS PRIOR TO OPERATION.

WARNING: THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM.

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A BRACE APPLICATION

Unfasten and loosen all straps and fasten hook/loop tabs back against the strap as shown. This will facilitate pulling the brace on.

While seated, pull on brace. The front hole must be positioned directly over the kneecap (patella), as shown. The plastic pressure plate should be positioned on the outside of the knee.

B TIGHTEN FOUR MAIN STRAPS

1st: lower rear; 2nd: upper rear;
3rd: lower front; 4th: upper front.

C TIGHTENING PATELLA BUTTRESS STRAPS

Bend the knee to 60°. While pushing down on the pressure plate, tighten the buttress straps around the kneecap. These straps should be snug.

NOTE: over-tightening will cause the medial (inside) hinge to pull forward.

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