POLAR CARE 300

CAST System

(Controlled Air Supply Tubing)

CAST System

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




5

The POLAR CARE 300 CAST System:

The POLAR CARE 300 CAST (Controlled Air Supply Tubing) System is designed to provide pneumatic power for the POLAR CARE 300 utilizing the standard oxygen or air supply commonly found in hospital rooms or in the Post Anesthesia Care Unit (PACU). The CAST System has a patented flow control and pressure relief valve (PRV) to provide the correct air flow and pressure for the POLAR CARE 300 to work optimally. The CAST System consists of tubing which connects to the oxygen or air supply, an in-line flow restrictor, a PRV, and tubing to connect the POLAR CARE 300. This system will enable the POLAR CARE 300 to function without any electrical supply and thereby eliminating the noise commonly heard from the standard electric air pump that is often in close proximity to the patient's head in the hospital setting. The CAST System also overcomes the problem of electrical outlets being installed upside down in the hospitals. The amount of compressed gas required to operate the POLAR CARE 300 is minimal.

Use instructions:

APPLY PAD ACCORDING TO INSTRUCTIONS INCLUDED WITH PAD.

- 1 PREPARE COOLER. Fill the cooler with ice, then with water to fill line, replace the lid.
- 2 CONNECT TO PAD. Attach the large blue hose from the top of the cooler to the pad. Push couplings together firmly until you hear a "click."
- 3 Place the cooler no more than 2 feet above or below the height of the pad.
- 4 Connect the 1/4" clear tubing of the CAST System to the compressed gas supply. Connect the blue air tube of the CAST System to the POLAR CARE 300.
- 5 Turn on the gas supply to at least 1 to 3 liters per minute. The CAST System will automatically regulate the gas flow and pressure to the correct settings for optimal function. The POLAR CARE 300 will cycle about every 18-25 seconds (a slightly audible release of gas will be heard from the top of the unit).

IF OXYGEN AND/OR AIR IS REQUIRED FOR OTHER PATIENT NEEDS, SIMPLY DISCONNECT THE CAST SYSTEM FROM THE COMPRESSED GAS SUPPLY.

To end session:

Disconnect the CAST System from the compressed gas supply. See step 4.

Replace ice and water before the ice melts completely.

Cleaning and maintenance:

After use, drain and dry pump with a soft cloth. Warm water and a mild detergent may be used occasionally to clean the cooler, lid and tubes.

Tips and problem solving:

If pad leaking is suspected, condensation is the probable cause. Some condensation on the lines, and pad is normal, particularly in humid environments. To help reduce this effect, provide additional insulating layers over the pad.

The couplings are self-sealing, however it is normal for a small amount of water to escape when disconnecting the couplings.

If the pad is not cooling properly, follow these steps:

- 1 Be sure the pad is not folded or kinked.
- 2 Be sure the pad's blue stem hose is not folded or kinked.
- 3 Make sure that the ice and water are at the proper level in the cooler.
- 4 Check air or gas flow. Be sure the air or gas supply is turned on to at least 1 liter per minute.
- 5 Check couplings for proper connection. Listen for "click" when connecting.

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