LPR™ KNEE BRACE

FITTING INSTRUCTIONS

ANLEGEANLEITUNG
ISTRUZIONI PER L'ADATTAMENTO
MISE EN PLACE
INSTRUCCIONES DE AJUSTE

BREG™

Breg, Inc.
2611 Commerce Way
Vista, CA 92081 U.S.A.
P: 800-321-0607
F: 800-329-2734
www.breg.com
Breg is an Orthofix company

AW-1.00351 Rev B 0311

	W A R N I N G S	
<p>WARNING: CAREFULLY READ FITTING INSTRUCTIONS AND WARNINGS PRIOR TO USE. TO ENSURE PROPER PERFORMANCE OF THE BRACE, FOLLOW ALL INSTRUCTIONS. FAILURE TO PROPERLY POSITION THE BRACE AND FASTEN EACH STRAP WILL COMPROMISE THE PERFORMANCE AND COMFORT.</p> <p>WARNING: FOR SINGLE PATIENT USE ONLY.</p> <p>WARNING: THIS PRODUCT IS NOT FOR USE ON SEDENTARY PATIENTS.</p> <p>WARNING: IF YOU EXPERIENCE INCREASED PAIN, SWELLING, SKIN IRRITATION, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.</p> <p>WARNING: THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM. CONSULT YOUR LICENSED HEALTH CARE PROFESSIONAL REGARDING SAFE AND APPROPRIATE ACTIVITY LEVEL WHILE WEARING THIS DEVICE.</p> <p>CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PRACTITIONER.</p> <p>CAUTION: CARE AND CLEANING OF THIS PRODUCT IS ESSENTIAL TO ITS CONTINUING STRENGTH AND PERFORMANCE. FOLLOW REMOVAL AND CLEANING PROCESS DESCRIBED IN INSTRUCTIONS.</p> <p>CAUTION: WHEN USING BRACES BILATERALLY, MODIFY ACTIVITY AS NECESSARY TO KEEP MEDIAL HINGES FROM CONTACTING EACH OTHER.</p> <p>CAUTION: WHEN USING AN LPR BRACE IT IS IMPORTANT TO KEEP DIRT AND GRIT OUT OF THE HINGE MECHANISM TO MAINTAIN PROPER HINGE MOBILITY AND BRACE FUNCTION.</p>		

Apply brace to leg (Example; right leg)

A. Unfasten all 5 straps, (brace is shown without straps in step A for clarity). Sit with the knee bent at a 30° angle. With the brace fully extended slide it onto the leg. Position the brace so that the center of the hinges align with the top of the patella.

Notice that the hinges may be positioned behind the bent knee, this is normal. As each strap is fastened, the hinges will move forward, until the brace is properly positioned on the leg with the center of the brace slightly back of the knee and midline of the leg (See Figure E).

B. With the brace on the leg, press the rigid shin cuff at the front of the brace flat against the leg and fasten strap 1 behind the calf. Next, fasten strap 2 in the back above the knee.

C. Be sure that the rigid shin cuff at the front of the brace is pressed flat against the leg and fasten strap 3, at the front of the brace.

If necessary, cut the strap to shorten it. Remove the numbered tab end, cut the strap to the desired length at the squared end. Replace the numbered tab end.

D. Next, fasten strap 4 in the back at the bottom of the brace; then fasten strap 5 in the back at the top of the brace. Be sure to fasten all of the straps in the proper order.


Tighten the straps snugly, but comfortably without over-tensioning any one strap. After 15 minutes of activity, sit with the knee bent at a 30° angle. Loosen and then retighten each strap, one at a time, in the proper order (this is necessary due to tissue swelling and compression). As you begin this process, be sure the hinge center aligns with the top of the patella (See Figure A). Improper strapping technique could result in discomfort. If you feel pressure on the thigh or shin, loosen the straps slightly. Press the rigid cuffs flat against the front of the thigh and shin, then retighten all the straps.

E. To remove the brace, unfasten all the straps on the back of the brace (1, 2, 4, and 5) and remove the brace off the front of the leg.

F. For PCL specific applications: (available on custom braces only). Press the rigid shin cuff at the front of the brace flat against the leg and fasten strap 1. Next fasten strap 2 in the back and above the knee making sure that the brace is slightly back of the knee and midline of the leg (see Figure F). Continue to follow the sequence as illustrated in Figure F until all six straps are applied, making sure to keep the center of the hinges at or slightly behind the midline.

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
Hinge range of motion adjustment

A. Brace is shipped with 10° hinge extension stops installed with additional flexion and extension stops in box. Stops are labeled in degrees:

Extension stops: 0°, 10°, 20°, 30° & 40°

Flexion stops: 45°, 60°, 75° & 90°

(brace is shipped with no flexion stops installed)

To change hinge stops, first remove the condyle pads held in place by hook/loop fasteners. Simply pull the pads off the hinge plates. To replace the pads be sure that they are positioned properly and firmly push them onto their hook/loop fasteners. Next, use a screwdriver to remove screws. Extension stops fit into the recessed slot on the face of the hinge plate, Flexion stops slide into the slot between the hinge plates from the back side. Remove existing stop, insert new stop, replace screw. Be sure that screws are seated and snug. Both medial (inside) hinge and lateral (outside) hinge stops must be identical in degrees.

Care and cleaning

Straps and pads must be removed and cleaned regularly depending upon your activity level and perspiration. Hand wash straps and frame pads in cold water with mild detergent, rinse, AIR DRY. Occasionally apply a dry lubricant to hinge especially after use in water.

B. The frame pads are held in place by hook/loop fasteners. Carefully pull the pads off the brace frame. To replace the pads be sure that they are positioned properly and firmly push them onto their hook/loop fasteners.

The knee hinge condyle pads are held in place by hook/loop fasteners. Simply pull the pads off the hinge plates. To replace the pads be sure that they are positioned properly and firmly push them onto their hook/loop fasteners.

C. If you require a more snug fit of the brace at the knee, shims can be installed behind the condyle pads. Pull the condyle pads off the hinge plates, then place the shim with the loop side against the inner hinge plate, then replace the pad over the shim. Be sure the pads are positioned properly and firmly push them onto their hook/loop fasteners.

D. To remove the brace straps, loosen the hook/loop tabs at each strap end and feed the straps through the D rings on the medial (inside) side of the brace until the strap is free. Be sure to replace the straps according to their strap number locations, See illustration.

E. Strap 3, (the front strap) must be replaced as shown. Both ends of the strap will thread through the D rings from the inside and fold over onto the outside. The strap end with the rectangular piece of hook folds across the front side of the strap first. Then the strap end with the plastic tab closes over the front of the strap.

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