# X2K<sup>®</sup> ADJUSTABLE HINGE, COUNTERFORCE PLUS and X2K-OA FITTING INSTRUCTIONS

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



**WARNINGS**

- CAREFULLY READ FITTING INSTRUCTIONS AND WARNINGS PRIOR TO USE.
- FOR SINGLE PATIENT USE ONLY.
- IF YOU EXPERIENCE INCREASED PAIN, SWELLING, SKIN IRRITATION, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.
- THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO ESSENTIAL PARTS OF A SAFE TREATMENT PROGRAM. CONSULT YOUR LICENSED HEALTH CARE PROFESSIONAL SHOULD YOU HAVE ANY QUESTIONS REGARDING THE SAFE APPLICATION OF THIS PRODUCT RELATED TO YOUR MEDICAL HISTORY.

CAUTION: CARE AND CLEANING OF THIS PRODUCT IS ESSENTIAL TO ITS CONTINUING STRENGTH AND PERFORMANCE. FOLLOW REMOVAL AND CLEANING PROCESS DESCRIBED BELOW.

CAUTION: WHEN USING BRACES BILATERALLY, MODIFY ACTIVITY AS NECESSARY TO KEEP MEDIAL HINGES FROM CONTACTING EACH OTHER.

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## Adjustable Hinge Instructions

### A ADJUSTMENT FOR LOOSE FITTING PADS

Use tool provided to turn upper and lower adjustment screws on sides of medial and lateral hinges clockwise. This will move condyle pads inward, tighter against knee. Adjust for comfortable, consistent fit.

### B ADJUSTMENT FOR TIGHT FITTING PADS

Turn upper and lower adjustment screws on sides of medial and lateral hinges counter clockwise. This will move condyle pads outward, away from knee.

### C ADJUSTMENT IF TOO LOOSE OR TIGHT AT UPPER OR LOWER EDGE OF PADS

Condyle pads may be tilted to fit contour of the knee. Upper or lower edge of pads may be adjusted to fit more tightly or loosely by turning upper or lower adjustment screws. Turn screws counter clockwise to move pad edge outward, away from knee. Turn screws clockwise to move pad edge inward, tighter against knee.

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
## CounterForce Plus and X2K-OA Unicompartmental Osteoarthritis Instructions

### A ADJUST THE LATERAL (OUTSIDE) CONDYLE PAD

With brace on, walk 5 minutes to position brace properly and ensure that lateral condyle pad is applying lateral force to knee. If more lateral force is required to relieve symptoms turn upper and lower adjustment screws on side of hinge clockwise. This will move condyle pad inward (toward knee), applying more lateral force to knee.

- 1 Upper and lower adjustment screws can be adjusted differently to accommodate various knee contours.
- 2 For more comfort, different adjustments of upper and lower screws will distribute lateral force more evenly.
- 3 Make sure large gel condyle pad is attached to the lateral (outside) condyle basket.

### B FIT THE MEDIAL (INSIDE) CONDYLE PAD

With brace properly fitted and lateral (outside) condyle pad adjusted for optimal comfort and pain relief, turn adjustment screws on side of medial hinge clockwise until condyle pad contacts knee with firm, comfortable pressure. The medial condyle pad can also be adjusted to accommodate various knee contours and maximize the pad's contact area.

**Important Note:** When wearing brace for an extended period you may experience discomfort caused by pressure from medial or lateral condyle pads. Back off adjustment screws to relieve pressure. Remove brace during periods of inactivity.

### Fitting instructions for lateral compartment symptoms

The procedure for lateral compartment symptoms is exactly the reverse of instructions for medial compartment symptoms (above). Ensure that brace fit and condyle pad adjustments are correct, readjust lateral condyle pad for firm, comfortable fit. Make sure large gel condyle pad is attached to the medial (inside) condyle basket.

### Fitting Instructions for gain or loss in leg girth:

#### C ADJUSTING CONDYLE PADS INWARD OR OUTWARD


If there is loss in girth or gain in girth (swelling), pressure against knee from condyle pads may be too great or insufficient. Condyle pads may be adjusted medially or laterally (inward or outward) to accommodate change in girth.

##### 1 ADJUSTMENT FOR LOSS IN GIRTH:

Use tool provided to turn upper and lower adjustment screws on sides of medial and lateral hinges clockwise. This will move condyle pads inward, tighter against knee.

##### 2 ADJUSTMENT FOR GAIN IN GIRTH (SWELLING):

Use tool provided to turn upper and lower adjustment screws on sides of medial and lateral hinges counter clockwise. This will move condyle pads outward, away from knee.

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