# Freestyle™ OA Brace

## FITTING INSTRUCTIONS

### for Medial OA application

ANLEITUNG ZUM ANLEGEN  
ISTRUZIONI PER L'ADATTAMENTO  
MISE EN PLACE  
INSTRUCCIONES DE AJUSTE

Breg, Inc.  
2885 Loker Ave. East  
Carlsbad, CA 92010 U.S.A.  
P: 800-321-0607  
F: 800-329-2734  
www.breg.com  
Breg is an Orthofix company

AW 1.93177 REV A 0911

EC REP

E/U authorized representative  
MDSS GmbH  
Schiffgraben 41  
D-30175 Hannover  
Germany

**⚠ WARNINGS ⚠**

**WARNING:** CAREFULLY READ FITTING INSTRUCTIONS AND WARNINGS PRIOR TO USE. TO ENSURE PROPER PERFORMANCE OF THE BRACE, FOLLOW ALL INSTRUCTIONS. FAILURE TO PROPERLY POSITION THE BRACE AND FASTEN EACH STRAP WILL COMPROMISE PERFORMANCE AND COMFORT.

**WARNING:** FOR SINGLE PATIENT USE ONLY.

**WARNING:** IF YOU EXPERIENCE INCREASED PAIN, SWELLING, SKIN IRRITATION, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.

**WARNING:** THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM. CONSULT YOUR LICENSED HEALTH CARE PROFESSIONAL REGARDING SAFE AND APPROPRIATE ACTIVITY LEVEL WHILE WEARING THIS DEVICE.

**CAUTION:** CARE AND CLEANING OF THIS PRODUCT IS ESSENTIAL TO ITS CONTINUING STRENGTH AND PERFORMANCE. FOLLOW REMOVAL AND CLEANING PROCESS DESCRIBED IN INSTRUCTIONS.

**CAUTION:** WHEN USING BRACES BILATERALLY, MODIFY ACTIVITY AS NECESSARY TO KEEP MEDIAL HINGES FROM CONTACTING EACH OTHER.

**CAUTION:** WHEN USING THE FREESTYLE OA BRACE IT IS IMPORTANT TO KEEP DIRT AND GRIT OUT OF THE HINGE MECHANISM TO MAINTAIN PROPER HINGE MOBILITY AND BRACE FUNCTION.

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


**PRECAUCIÓN:**

- A** Prior to application, set the hinge to 0.
- B** Sit with the knee bent at 30°
- C** Position the brace by placing the condyle pad on the outside of the knee.
- D** Maintain the brace position with one hand while using the other hand to pull tab #1 across the thigh above the knee and fasten securely.
- E** Pull tab #2 across the shin below the knee and fasten securely.

- F** Pull tab #3 across the lower shin and fasten securely.
- G** Pull tab #4 across the upper thigh and fasten securely.
- H** Pass strap #5 through the plastic ring on the outside of the thigh.
- I** To avoid rotation of the brace during application, secure the brace with one hand, as shown, while tensioning strap #5. Lay the strap onto itself and fasten securely. Repeat for strap #6.
- J** Once the brace is properly applied, pain relief may be achieved by turning the thumbwheel hinge dial clockwise to increase pressure to the outside of the knee.

Scan QR Code for  
Instruction Video


 [breg.com/foa](http://breg.com/foa)

**A**

**Instruction Video**

To view a video demonstrating how to use the brace, scan the QR code or go to the website: [www.breg.com/foa](http://www.breg.com/foa)


**B**

**Hinge Adjustment**

With the brace properly applied, OA pain symptoms can be relieved by adjusting the thumbwheel hinge dial. Turn the thumbwheel dial clockwise to increase pressure and counterclockwise to reduce the pressure on the outside of the knee. The hinge settings window displays settings ranging from 0 to 10mm. The display setting will increase or decrease as you adjust the thumbwheel dial.

Make sure a large gel condyle pad is attached inside the brace.


Always follow your medical professional's recommendation regarding the proper hinge setting.

**C**

**Product diagram**

- |                    |                                  |
|--------------------|----------------------------------|
| <b>1:</b> Tab #1   | <b>7:</b> Condyle pad            |
| <b>2:</b> Tab #2   | <b>8:</b> Finger Pockets         |
| <b>3:</b> Tab #3   | <b>9:</b> Hinge Settings Window  |
| <b>4:</b> Tab #4   | <b>10:</b> Hinge Thumbwheel Dial |
| <b>5:</b> Strap #5 | <b>11:</b> Plastic Rings         |
| <b>6:</b> Strap #6 |                                  |

**D**

**Condyle Pad Adjustment**

If you require a more snug fit of the brace at the knee, shims can be installed to the condyle pad.

1. Remove the condyle pad cover and carefully peel off the condyle pad that is attached with hook/loop.
2. Attach the desired number of shim(s) to the condyle pad.
3. Re-attach the condyle pad to the brace.
4. Turn the Hinge Thumbwheel Dial clockwise until the desired pressure is reached.

Important: Make sure the condyle pad is properly aligned with the hinge thumbwheel dial before applying the brace to the knee.

**Care and Cleaning**

Gently hand wash the material in cold water with a mild detergent, rinse, lay flat and air dry. Apply a dry lubricant such as a Teflon spray to the hinge after washing. **DO NOT** machine wash or tumble dry.

**Fitting Tips**

- A properly applied brace should fit snug on all areas the knee.
- Use the finger pockets to help fasten tabs #1 – #4.
- Avoid rotation of the brace during application of straps #5 and #6. Hold the brace securely with one hand at the hinge near the plastic ring while tensioning the straps with the other hand.