BOA Classic Spine Brace

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BOA Chairback Instructions on reverse side

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AW-1.90110 REV 8 0709

WARNING: CAREFULLY READ USE/CARE INSTRUCTIONS AND WARNINGS PRIOR TO USE.

WARNING: THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM.

WARNING: IF YOU EXPERIENCE INCREASED PAIN, SWELLING, SENSATION CHANGES, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.

WARNING: IF YOU EXPERIENCE ANY SKIN IRRITATION OR ADVERSE SKIN REACTION WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR HEALTH CARE PROFESSIONAL.

WARNING: IF ANY FRAYING OR DETERIORATION OF THE CORDING IS NOTICED, DISCONTINUE USE.

CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH PRACTITIONER.




CAUTION: FOR SINGLE PATIENT USE ONLY.

- Position and fasten the pull-handle as close as possible to the compression panels, then fully extend the cording between the 2 halves of the brace.
- IF USING THE LORDOTIC FOAM INSERTS: Slide the appropriate insert (7°, 10°, or 12°) into the elastic pockets located on the inside of the brace halves. Adjust for comfort.
- Position the anterior panel over the center of the patient's abdomen with bottom of the panel just above the pubic bone.
NOTE: The patient will need assistance holding the anterior panel in position in order to complete steps 4 & 5!
- Position the brace behind the patient's back with the interior lining facing the body. For proper fit and compression, the opening in the middle of the brace is centered over the spine. The pull-handle will be on the patient's left side.
- Wrap the left side (with pull-handle) of the brace around the waist and secure it to the anterior panel by attaching the hook strip on the panel to the corresponding loop strip on the inside of the brace.
- Next, wrap the right side around the waist and fasten.
- To apply compression, pull the (pull-handle) away from the body until the desired level of compression is achieved, fasten the pull-handle to the contact strip.

CARE INSTRUCTIONS

Hand wash in cold water using mild soap. Rinse thoroughly, dry with towel then air dry completely. If not rinsed thoroughly, soap residue may cause skin irritation and material breakdown.

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


A. Position and fasten the pull-handle as close as possible to the compression panels, then fully extend the cording between the 2 halves of the brace.

B. Place the chairback on a flat surface with the foam interior facing down. Open the straps on each side of the chairback, place the brace over the chairback so that the middle of the brace is centered, re-fasten the straps to secure the brace to the frame.

C. Position the anterior panel over the center of the patient's abdomen with bottom of the panel just above the pubic bone.
NOTE: The patient will need assistance holding the anterior panel in position in order to complete steps D & E !

D. Position the brace, with chairback attached, behind the patient's back. For proper fit and compression, the opening in the middle of the brace is centered over the spine. The pull-handle will be on the patient's left side.

E. Wrap the left side (with pull-handle) of the brace around the waist and secure it to the anterior panel by attaching the hook strip on the panel to the corresponding loop strip on the inside of the brace.

F. Next, wrap the right side around the waist and fasten.

G. To apply compression, pull the (pull-handle) away from the body until the desired level of compression is achieved, fasten the pull-handle to the contact strip.

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