Svenska

1. ÖPPNA SPÄNNENA

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INDICATIONS: For increased support for patients who have continued symptoms of giving way, poor quadriceps or hamstring strength, or a desire for an early return to sports activities. Indicated for increased support of knees following injury to or reconstruction of the anterior or posterior cruciate ligament, medial or lateral collateral ligament, or menisci either alone or in combinations when properly adjusted.

CONTRAINDICATIONS:

Contraindicated in patients with triple varus knees that have damage to the postero-lateral structures combined with stretched lateral structures and varus thrust.

WARNING:

CAUTION: For Single Patient Use Only.
CAUTION: Activity modification and proper rehabilitation are an essential part of a safe treatment program. If any additional pain or other symptoms occur while using this device, discontinue use and seek medical attention.

CAUTION:

Do not leave this device in a car on a hot day, or in any other place where the temperature may exceed 150°F (65° C.) for any extended period of time.

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


English

1. UNFASTEN BUCKLES

Release all four straps on the back side of the brace by squeezing the gray levers on either side of the buckle and withdrawing the over-center buckle. Loosen the front straps.

2. POSITION THE BRACE

With the knee bent at 45°, lay the brace on the leg and align the hinge as shown above. Position the center of the hinge on an imaginary line that goes across the top of the knee cap. Push the hinges posteriorly until the hinges are in the center of the knee joint when viewed from the side **(A)**.

3. ANCHOR THE BRACE

Wrap the upper rear calf strap around the calf as high as possible over the calf muscle. Insert the over-center tab into the base. Push the tab until you hear the hinge click on both sides **(B)**.

Adjust the strap by releasing the tab and pull the strap until a snug and comfortable fit is achieved. All straps have additional length adjustment at the opposite end of each strap. See step 5 **(C)**.

(Z-12-D version: Feed the tab end of the strap through the D-ring. Tension the strap until a snug and comfortable fit is achieved.)

4. FASTEN LOWER CALF STRAP

Insert the over-center tab into the buckle base on the lower calf shell of the brace. Push the tab in until both sides click locked. Adjust strap until a comfortable fit is achieved. (If the strap tab

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5. ADJUST STRAP LENGTH

All straps have additional length adjustment at the end opposite the tab. Release the strap and adjust the length of the strap using the hook and loop fastener **(D)**.

6. FASTEN REAR THIGH STRAPS

Following the same instructions as in Step 4, fasten the lower and then upper rear thigh strap.

7. ADJUST FRONT TIBIA STRAP

Tighten the front tibia strap. If the strap overlaps too far, see Step 6 for adjustment instructions. Center the strap pad **(E)**.

8. ADJUST FRONT THIGH STRAP (PCL/COMBINED INSTABILITY VERSIONS)

Tighten the front thigh strap. (Note: Over tightening the front thigh strap can cause the brace not to function properly.) Stand up and walk or jog around to see how your new Z-12 brace feels. Straps should not be so tight that blood circulation is restricted. Re-adjust as needed. If this is your first time to wear a functional knee brace, it might take a few hours of wear to get used to how it feels **(F)**.

9. ADJUST OA CONDYLE PAD

Press down on gray tab and slide into position for desired unloading. Adjust for comfort **(G)**.


Patient Information

SPECIAL ADJUSTMENT FOR ANTERIOR CRUCIATE LIGAMENT INJURIES

In order to provide increased resistance against tibial subluxation (forward movement of the lower leg relative to the upper leg) when the leg is extended or to reduce force on a reconstructed ACL, it may be desirable to pretension the straps to provide a “shear force” across the knee by increasing the tension on the front calf strap on the tibia and the lower rear thigh strap on the back of the thigh. All straps for ACL patients should be adjusted with the leg bent at 60°-80°. Loosen the front thigh strap and increase the tension on the lower rear thigh strap and, if necessary, the front shin strap to provide a pre-loaded strap force that helps prevent anterior tibial subluxation.

SPECIAL ADJUSTMENT FOR POSTERIOR CRUCIATE LIGAMENT INJURIES

In order to provide increased resistance against posterior tibial subluxation (backward movement of the lower leg relative to the upper leg) as the leg bends or to reduce force on a PCL reconstruction, it may be desirable to pretension the straps to provide a “shear force” across the knee by increasing the tension on the upper rear calf strap and the front thigh strap. All straps for PCL patients should be adjusted with the leg in full extension. Loosen the lower rear thigh strap and the front calf strap and increase the tension on the front thigh strap and the upper rear calf strap to provide a pre-loaded shear force that prevents posterior tibial subluxation. NOTE: If needed, apply the included spacers behind the condyle pads to achieve proper condular pressure on the medial and/or lateral side of the knee.

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© 2015 Breg, Inc. All rights reserved.