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WARNING: Federal Law Restricts this Device to Sale by or on the Order of a Licensed Health Care Practitioner.

WARNING: For Single Patient Use Only.

CAUTION: Activity modification and proper rehabilitation are an essential part of a safe treatment program. If any additional pain or other symptoms occur while using this device, discontinue use and seek medical attention.

CAUTION: This device will not prevent or reduce all injuries.

CAUTION: Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions. Failure to properly position the brace and fasten each strap will compromise performance and comfort.

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CAUTION: This device has been developed for the indications listed above, and may be damaged by unusual activities such as sports, falls, or other accidents, or exposure of the hinges to debris or liquids. Do not leave this device in a car on a hot day, or in any other place where the temperature may exceed 150°F (65° C) for any extended period of time. To avoid the risk of accident and possible injury to the patient, discontinue use if such an incident has occurred.

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







ENGLISH

- Unfasten Buckles**
Release all four straps on the back side of the brace by squeezing the gray levers on either side of the buckle and withdrawing the over-center buckle. Loosen the front strap(s).
- Positioning the Brace**
ACL: With the knee bent at 60°, locate the top of the knee cap.
PCL: With the knee bent at 10°, locate the top of the knee cap.
Position the center of the hinge on an imaginary line that goes across the top of the knee cap. Push the hinges posteriorly until the hinges are in the center of the knee joint when viewed from the side. (A)
- Anchoring the Brace**
Wrap the upper rear calf strap around the calf as high as possible over the calf muscle. (B) AXIOM-D: Place the tab end of the strap through the D-ring. Pull the tab until a snug and comfortable fit is achieved.
AXIOM: Insert the over-center buckle into the base. Push the buckle until you hear the gray levers click on both sides.
Adjust the strap by releasing the tab and pull the strap until a snug and comfortable fit is achieved. (C)
- Fastening Lower Calf Strap**
Insert the over-center buckle into the base on the lower calf shell of the brace. Push the buckle in until you hear the gray levers on click on both sides. Adjust strap until a comfortable fit is achieved. (If the strap tab reaches its limit without achieving the desired fit, see strap length adjustment in figure (E).)
- Fastening Rear Thigh Straps**
Following the same instructions as in step 4, fasten the lower and then upper rear thigh strap.
- Adjusting Front Tibia Strap**
Tighten the front tibia strap. If the strap overlaps too far, see length adjustment instructions. (D)
- Adjust Strap Length**
All straps have additional length adjustment at the end opposite the tab. Release the strap and adjust the length of the strap using the hook and pile fastener. Center the strap pad. (E)
- Adjust Front Thigh Strap (PCL/Combined Instability Versions)**
Tighten the front thigh strap. **Note:** Over tightening the front thigh strap can cause the brace not to function properly. Stand up and walk or jog around to see how your new AXIOM or AXIOM-D brace feels. Straps should not be so tight that blood circulation is restricted. Re-adjust as needed. If this is your first time to wear a functional knee brace, it might take a few hours of wear to get used to how it feels. (F)
- Adjusting OA Condyle Pad**
While sitting or with the brace off the leg, press down on the gray tab on dial and slide into position for desired unloading. Adjust for comfort. (G)

INDICATIONS: For increased support for patients who have continued symptoms of giving way, poor quadriceps or hamstring strength, or a desire for an early return to sports activities.
-AXIOM: Indicated for increased support of knees following injury to or reconstruction of the anterior or posterior cruciate ligament, medial or lateral collateral ligament, or menisci either alone or in combinations when properly adjusted.
-AXIOM-D: Indicated for increased support of knees following injury to or reconstruction of the anterior or posterior cruciate ligament when properly adjusted.
-AXIOM/AXIOM-D ADJ OA: Indicated for patients with combined instability as well as mild to moderate osteoarthritis. Not recommended for patients with abnormal amount of redundant tissue.
CONTRAINDICATIONS: Contraindicated in patients with triple varus knees that have damage to the postero-lateral structures combined with stretched lateral structures and varus thrust.

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© 2015 Breg, Inc. All Rights Reserved