## Swedish

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
### NOTICE TO USERS

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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**WARNING:** For Single Patient Use Only.

**CAUTION:** United States Federal Law Restricts this Device to Sale by or on the Order of a Licensed Health Care Practitioner.

**CAUTION:** Activity modification and proper rehabilitation are an essential part of a safe treatment program. If any additional pain or other symptoms occur while using this device, discontinue use and seek medical attention.

**CAUTION:** This device will not prevent or reduce all injuries.

**CAUTION:** Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions. Failure to

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
## English

### 1. UNFASTEN BUCKLES

Release all four straps on the back side of the brace by squeezing the gray levers on either side of the buckle and withdrawing the over-center buckle. Loosen the front strap(s).

### 2. POSITIONING THE BRACE

**ACL:** With the knee bent at 60°, locate the top of the knee cap.

**PCL:** With the knee bent at 10°, locate the top of the knee cap.

Position the center of the hinge on an imaginary line that goes across the top of the knee cap. Push the hinges posteriorly until the hinges are in the center of the knee joint when viewed from the side. **(A)**

### 3. ANCHORING THE BRACE

Wrap the upper rear calf strap around the calf as high as possible over the calf muscle. **(B)**

**AXIOM-D ELITE:** Place the tab end of the strap through the D-ring. Pull the tab until a snug and comfortable fit is achieved.

**AXIOM ELITE:** Insert the over-center buckle into the base. Push the buckle until you hear the gray levers click on both sides.

Adjust the strap by releasing the tab and pull the strap until a snug and comfortable fit is achieved. **(C)**

### 4. FASTENING LOWER CALF STRAP

Insert the over-center buckle into the base on the lower calf shell of the brace. Push the buckle in until you hear the gray levers click on both sides. Adjust strap until a comfortable fit is achieved. (If the strap tab reaches its limit without achieving the desired fit, see strap length adjustment in figure **(E)**.)

**INDICATIONS:** For increased support for patients who have continued symptoms of giving way, poor quadriceps or hamstring strength, or a desire for an early return to sports activities.

**AXIOM ELITE:** Indicated for increased support of knees following injury to or reconstruction of the anterior or posterior cruciate ligament, medial or lateral collateral ligament, or menisci either alone or in combinations when properly adjusted.

**AXIOM-D ELITE:** Indicated for increased support of knees following injury to or reconstruction of the anterior or posterior cruciate ligament when properly adjusted.

### 5. FASTENING REAR THIGH STRAPS

Following the same instructions as in step 4, fasten the lower and then upper rear thigh strap.

### 6. ADJUSTING FRONT TIBIA STRAP

Tighten the front tibia strap. If the strap overlaps too far, see length adjustment instructions. **(D)**

### 7. ADJUST STRAP LENGTH

All straps have additional length adjustment at the end opposite the tab. Release the strap and adjust the length of the strap using the hook and pile fastener. Center the strap pad. **(E)**

**8. ADJUST FRONT THIGH STRAP (PCL/ COMBINED INSTABILITY VERSIONS)** Tighten the front thigh strap. Note: Over tightening the front thigh strap can cause the brace not to function properly. Stand up and walk or jog around to see how your new **AXIOM ELITE** or **AXIOM-D ELITE** brace feels. Straps should not be so tight that blood circulation is restricted.

Re-adjust as needed. If this is your first time to wear a functional knee brace, it might take a few hours of wear to get used to how it feels. **(F)**

Re-adjust as needed. If this is your first time to wear a functional knee brace, it might take a few hours of wear to get used to how it feels. **(F)**

### 9. CLEANING

Pad and straps may be hand washed in cold water and mild detergent. Rinse and lay flat to air dry. Do not machine wash or tumble dry.

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