English

1. Apply the Full Boot Cuff
(For boots with a calf cuff (open toe), proceed to step 2).

a. Slip the foot into the bootie with the toes beneath the tongue. Slide the heel of the foot all the way to the back of the bootie.

b. Close the foam on the forefoot.

2. Secure the Cuff
a. Wrap the foam section around the leg and tuck in the long end of the foam cuff.

b. Adjust the tension on the wrap, and close the hook and loop closure.

c. If the cuff is too long, trim with scissors.

3. Apply the Boot




Make sure the plastic covers are still on the boot uprights. Slip the foot into the boot, position the foot with the ankle bones centered between the uprights.

4. Secure the Uprights

a. While bearing weight and with the foot pointed directly forward, align the uprights with the sides of the leg parallel to the leg bones. When the uprights are properly aligned with the leg, remove the plastic from the uprights and press the uprights against the sides of the leg to engage the hook material to the bootie or cuff.

b. Save the plastic covers for reapplication after washing the bootie or cuff.

5. Secure the Forefoot Straps
While standing, wrap the forefoot straps over the top of the foot, through the D-rings and pull back to tension. If necessary, trim the excess strap length with scissors.

Deutsch

6. Adjust the Rear Plate
(For boots without a rear plate, proceed to step 8)

If necessary, reapply the plastic plate containing the straps so that the top of the plate is positioned about an inch below the top edge of the foam cuff. Make sure the plate is centered (left to right) on the cuff at the back of the calf. The straps should remain attached to the plate. NOTE: If necessary, the plate may be repositioned higher or lower to allow the straps to align with the attachments on the uprights.

7. Secure the Leg Straps

a. While still standing, beginning at the ankle wrap both strap ends firmly forward to engage the straps with the hook material on the outside of the uprights.

b. Loop the end of the strap through the D-ring and pull back to tension. Press hook end to strap to engage.

c. Repeat steps A and B for the remaining straps.

NOTE: For mid-calf versions there is only one remaining leg strap.

...Instructions Continue on Reverse Side

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