

## HORIZON™ 456 TLSO

Distributed by:

Breg, Inc.  
2885 Loker Ave East  
Carlsbad, CA 92010 USA  
(800) 321-0607  
www.breg.com

### ENGLISH

#### SIZE ADJUSTMENT

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
### ENGLISH

#### PREPARATION

Place the pull-tabs on the belt 2 inches (5 cm) from the plastic in the back of the Horizon™ 456. Pull the ends of the belt away from each other ensuring that the belt is fully extended.

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
**ENGLISH**  
**APPLICATION**

a. Apply the fully extended brace making sure the back panel is centered and low on the back. Apply the brace tightly around the abdomen, the right side overlapping the left.

Note: When applied properly, both the upper and lower edge of the brace should be snug. If there is any gap, reapply the brace and angle the ends to insure a snug fit at both the upper and lower edge. For certain patients, the ends must be steeply angled.

b. Slip each thumb into the pull tabs. Extend the tabs out and away from the brace. When desired compression is achieved, place the hook side of the pull tabs down onto the loop material to secure.

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
**ENGLISH**  
**ADJUSTMENT**

a. Adjust the vertical height with the posterior adjustment panel. The yoke should sit at level C7 on patient.

b. Pull and tighten the shoulder straps until properly fit.

c. Fasten and tighten the chest buckle to draw the shoulder straps away from the armpits which may be more desirable for some patients.

Note: After initial fitting, have patient sit, then stand, and re-adjust fitting as necessary.

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**LIMITED WARRANTY**  
Aspen Medical Products, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products in the event of breach of warranty shall be to repair or replace the defective product or part(s).  
Aspen Medical Products shall have no obligation under this limited warranty in the event:  
(a) The product was not purchased from Aspen Medical Products or through its authorized channels of distribution;  
(b) The product is altered;  
(c) Any parts not supplied by Aspen Medical Products are inserted into the product; or  
(d) The product is not used in accordance with the Aspen Medical Products, Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OF, OR DEFECT IN, THE PRODUCT. ASPEN MEDICAL PRODUCTS SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS.  
This warranty gives you specific legal rights and you may have other rights which vary from state to state.  
Warning: To ensure proper use, please review all material.  
A listing of authorized distributors in your area is available upon request.

**READ INSTRUCTIONS BEFORE USE.** Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations.

THIS PRODUCT IS INTENDED FOR APPLICATION BY LICENSED HEALTH CARE PRACTITIONERS AS DIRECTED BY A PHYSICIAN OR OTHER QUALIFIED MEDICAL AUTHORITY. PERMANENT INJURY MAY RESULT FROM USE OTHER THAN AS DIRECTED BY A PHYSICIAN OR OTHER MEDICAL AUTHORITY. ASPEN MEDICAL PRODUCTS SPECIFICALLY DISCLAIMS ANY OBLIGATION OR LIABILITY FOR INJURY OR DAMAGE BY REASON OF ANY APPLICATION OF THIS PRODUCT. PLEASE SEE THE ACCOMPANYING LIMITED WARRANTY FOR FURTHER INFORMATION IN THIS REGARD.