
BREG

LEGACY THRUSTER

INSTRUCTIONS FOR USE

CP020256 Rev E 09/15

CAUTION: Activity modification and proper rehabilitation are an essential part of a safe treatment program. If any additional pain or other symptoms occur while using this device, discontinue use and seek medical attention.

CAUTION: This device is not intended to prevent injury, or reduce or eliminate the risk of re-injury.

CAUTION: This device has been developed for the indications listed above, and may be damaged by unusual activities such as sports, falls, or other accidents, or exposure of the hinges to debris or liquids. Do not leave this device in a car on a hot day, or in any other place where the temperature may exceed 150F (65 C.) for any extended period of time. To avoid the risk of accident and possible injury to the patient, discontinue use if such an incident has occurred.

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This device is offered for sale by or on the order of a physician or other qualified medical professional. This device is not intended for reuse on a second patient. This device is not intended to prevent injuries or reduce or eliminate the risk of re-injury.

INDICATIONS: The Legacy Thruster brace can be configured to address medial or lateral osteoarthritis. The brace is designed to decrease joint pain by creating a space between the bones on the affected side of the knee.

CONTRAINDICATIONS: The Legacy Thruster is contraindicated to treat medial osteoarthritis with a lateral osteoarthritis brace or to treat lateral osteoarthritis with a medial osteoarthritis brace.

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© 2015 Breg, Inc. All Rights Reserved
2885 Loker Ave E
Carlsbad, CA 92010 USA

Toll Free Tel: 1-800-321-0607

Local: +1-760-795-5440

Fax: +1-760-795-5295

www.Breg.com

CE

EC REP

EU authorized representative

MDSS GmbH

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





D-30175 Hannover

Germany

LEGACY THRUSTER

ENGLISH

1 Unfasten Buckle

Unfasten all of the buckle ends by pushing the gray buttons on each side of the buckle **(A)**. If this is the first time the brace is applied, the brace should be at its lowest setting **(B)**.

2 Apply the Brace

Apply the upper calf strap by wrapping the strap behind the calf muscle and around to the front of the shell to engage the buckle **(D)**. Adjust the length of the strap using the hook and pile fastener at the buckle end of the strap.
MEDIAL COMPARTMENT OA: Position the hinge on the **INSIDE** of the leg and in line with the top of the knee cap.
LATERAL COMPARTMENT OA: Position the hinge on the **OUTSIDE** of the leg and in line with the top of the knee cap.

3 Fasten Upper Calf Strap

Apply the upper calf strap by wrapping the strap behind the calf muscle and around to the front of the shell to engage the buckle **(D)**. Adjust the length of the strap using the hook and pile fastener at the buckle end of the strap.

4 Fasten Lower Calf Strap

Apply the lower calf strap in the same manner **(E)**.

5 Fasten Lower Thigh Strap

Wrap the lower thigh strap from the front of the thigh around the back of the leg and engage buckle on the opposite side of the shell. Adjust the length of the strap using the fastener at the buckle end of the strap.

6 Fasten Upper Thigh Strap

Wrap the upper thigh strap around the back of the leg to engage the buckle. Adjust the strap using the fastener at the buckle end of the strap.

7 Hinge Adjustment

Stand and walk in the brace to see how it feels on the leg. Some may experience pain relief prior to adjusting the hinge, due to varying leg angulations. The brace only applies force as it approaches full extension. It can only be adjusted while seated with the leg bent at 90° in small increments. If more force is required, sit down and adjust the dial to an increased force setting. Walk again to test for pain relief **(F)**.

NOTE: Only adjust the hinge enough for pain relief. Walk after each adjustment. Too much force may cause discomfort and excessive muscle soreness. Once you reach pain relief, don't adjust the hinge any further. It takes at least a week for the leg to get used to having pressure applied to it. It is important that the straps closest to the knee are snug, or you will lose brace effectiveness. If there is any pinching from the straps, bend the leg to 90° and readjust the straps for comfort.

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