# Patella Tracking Orthosis (PTO) Knee Brace

## PTO-KNIESCHIENE

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
© 2023 Breg, Inc. All rights reserved.  
AirMesh is a registered trademark of Breg, Inc.

E/U authorized  
representative  
MDSS GmbH  
Schiffgraben 41  
D-30175 Hannover  
Germany

AW-1.08770 Rev J 10/23

## WARNING

**CAUTION:** FOR SINGLE PATIENT USE ONLY.

**CAUTION:** IF YOU EXPERIENCE INCREASED PAIN, SWELLING, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.

**WARNING:** CAREFULLY READ USE INSTRUCTIONS AND WARNINGS PRIOR TO USE.

**WARNING:** THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM.

## A BRACE APPLICATION

Unfasten and loosen all straps and fasten hook/loop tabs back against the strap as shown. This will facilitate pulling on the brace.

While seated, pull on brace. The front hole must be positioned directly over the kneecap (patella), as shown. The plastic pressure plate should be positioned on the outside of the knee.

## B TIGHTEN FOUR MAIN STRAPS

1st: lower rear; 2nd: upper rear;  
3rd: lower front; 4th: upper front.

## C TIGHTENING PATELLA BUTTRESS STRAPS

Bend the knee to 60°. While pushing down on the pressure plate, tighten the buttress straps around the kneecap. These straps should be snug.

NOTE: over-tightening will cause the medial (inside) hinge to pull forward.

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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