







Size	Waist Size
1	24" – 34" (61.0–86.4cm)
2	33" – 40" (83.8–101.6cm)
3	39" – 46" (99.1–116.8cm)
4	45" – 50" (114.3–127cm)
EXT.	50" + (127cm+)

ENGLISH

INTENDED USE/INDICATIONS:

The intended use of these lumbar products is to provide stabilization and motion restriction for the treatment of trauma related injury, acute post-surgical recovery, and the management of chronic conditions of the spine. Indicated for acute and chronic lower back pain, sprains/strains, lumbar disc displacement, deconditioned trunk musculature, facet syndrome, post-laminectomy syndrome, spinal stenosis, spondylolisthesis, spondylitis, and radiculopathy.

INTENDED USERS:

The intended user should be a licensed professional, the patient or patient's caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:

Not Applicable.

WARNINGS AND PRECAUTIONS:

- Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions.
- If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your licensed health care professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.
- This device is offered for sale by or on the order of a licensed health care professional.
- Federal law restricts this device to sale by or on the order of a licensed health care practitioner.
- The frequency and duration of use shall be prescribed by a healthcare professional.
- Follow application and care process described in instructions.
- For single patient use only.
- Do not use over open wounds.
- Do not use this device if damaged.

INSTRUCTIONS FOR USE: MEASUREMENT

1. Measure patient's waist and refer to sizing chart for corresponding size.

ADJUSTING THE WAIST BELT SIZE

Note: The initial configuration is set for size 2.

1. Open the brace and lay it face down.
2. Flip the belt to find sizing numbers.
3. Disengage hook and loop fastener and move each end to the desired size.
4. Repeat on both sides.

Note: Before trimming any material flip the belt over and wrap around patient's waist to estimate fit, then adjust accordingly.

APPLYING THE BRACE

Ensure pull handles are near the back of the nylon cover

5. Align pad centered with the spine and bottom edge over the lumbosacral joint.
6. Fasten brace by pulling right belt over left.
7. Adjust tension by pulling handles outward.
8. Attach handles to belt.

REMOVING THE BRACE

9. Detach pull handles and place near the back.
10. Unfasten front hook and loop to remove.

USE AND CARE:

Hand wash in 30°C (86°F) water with mild soap and rinse thoroughly. Lay flat to air dry completely before wearing. Note: if not rinsed thoroughly, residual soap can cause skin irritation or deterioration of the material.

MATERIAL CONTENT:

Airmesh®, Cotton, High Density Polyethylene (HDPE), Nylon, Polyester

WARRANTY:

Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:

Dispose of product according to local regulations.

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NOTICE TO USERS:

If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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