ENGLISH

INTENDED USE/INDICATIONS:

The intended use of these Lumbar products is to provide stabilization and motion restriction for the treatment of trauma related injury, acute post-surgical recovery, and the management of chronic conditions of the spine.

Indicated for bulging/herniated disc, decompressive procedures, deconditioned trunk musculature, degenerative disc disease, post-operative support, spinal stenosis, sprains/strains, spondylolisthesis, Sacroiliitis, Radiculopathy

INTENDED USERS:

The intended user should be a licensed professional, the patient, or patient's caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:

Not Applicable.

⚠️ WARNINGS AND PRECAUTIONS:

- Carefully read use instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions.
- If you experience increased pain, swelling, or any adverse reactions while using this product, immediately consult your medical professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.
- Follow application and care process described in instructions.
- For single patient use only.

INSTRUCTIONS FOR USE:

1. Adjust and secure the brace by moving both panels in or out to the desired size. The size markings are clearly labeled on each panel. Align and firmly secure each panel to hook and loop connector (A). If patient requires a larger size than Universal, a 4XL extension is available.
2. Position the brace so that it is centered and low on the back.
3. Pull belt panel ends forward, around the waist placing right belt panel on top of left belt panel to engage hook and loop.
4. Have patient slip their thumbs through the adjustment pull rings. Pull forward and around for secured adjustment.
5. Apply comfortably, and then secure pull rings by pressing firmly onto front of the brace. Adjust as needed.

USE AND CARE:

Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or tumble dry.

MATERIAL CONTENT:

Nylon, Polyester, Polypropylene, Polyurethane

WARRANTY:

Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:

Dispose of product according to local regulations.

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© 2025 Breg, Inc. All rights reserved.