




Size	Waist Size
1	26" – 32" (66–81.3cm)
2	31" – 37" (78.7–94cm)
3	36" – 43" (91.4–109.2cm)
4	42" – 49" (106.7–124.5cm)
5	48" – 55" (121.9–139.7cm)
6	54" – 60" (137.2–152.4cm)
EXT.	60" + (152.4cm+)


1


2


3


4

5

6

7

8

9

10


11


12

13

14

15

16

17

18

ENGLISH

INTENDED USE/INDICATIONS:

The intended use of the Pinnacle lumbar orthoses is to provide stabilization and motion restriction for the treatment of trauma related injury, acute post-surgical recovery, and the management of chronic conditions of the spine. Indicated for chronic low back pain, degenerative disc disease, severe DJD, radiculopathy, spondylosis, spinal stenosis, post laminectomy syndrome, spondylolisthesis, decompressive procedures, deconditioned trunk musculature, sprains/strains, bulging/herniated disc, facet syndrome, IDET procedures, SI dysfunction, sacroiliitis, post-operative support and fracture management.

INTENDED USERS:

The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:

Not Applicable.

⚠ WARNINGS AND PRECAUTIONS:

- Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions.
- If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your licensed health care professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.
- This device is offered for sale by or on the order of a licensed health care professional.
- For single patient use only.
- Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
- Follow the licensing healthcare practitioner's protocol regarding method, frequency and duration of administration.

INSTRUCTIONS FOR USE:

MEASUREMENT

1. Measure patient's waist and refer to sizing chart (above) for corresponding size.

ADJUSTING THE WAIST BELT

For Sizes 1–3:

2. Open brace and lay it face down.
3. Lift mesh cover and flip pulley plate.
4. Detach belt from pulley plate.
5. Slide belt through slot and fold to match size number.
6. Reattach belt and close mesh.

For Sizes 4–6:

7. Remove hook piece from pulley plate.
8. Move belt from inner to outer slot.
9. Slide and fold belt to match size number.
10. Reattach and secure pulley plate and mesh.

Repeat steps for both sides.

11. Re-attach the pull handles to the belt, then pull the belt ends away from each other to full extend the belt.

LORDOTIC PADS

12. Located inside the brace.

Can be removed or repositioned through the center hole in the back pad.

ADJUSTING THE BRACE

13. Align back pad with spine and lumbosacral joint.
14. Wrap right belt over left and fasten.
15. Pull handles to adjust tension.
16. Attach handles to belt.

REMOVING THE BRACE

17. Detach handles and place near back.
18. Undo front hook-and-loop closure.

USE AND CARE:

Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or tumble dry.

MATERIAL CONTENT:

Airmesh®, Cotton, High Density Polyethylene (HDPE), Neoprene, Nylon

WARRANTY:

Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:

Dispose of product according to local regulations.

NOTICE TO USERS: If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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