




Size	Waist Size
1	22" – 28" (55.9–71.1cm)
2	27" – 33" (68.6–83.8cm)
3	32" – 39" (81.3–99.1cm)
4	38" – 45" (96.5–114.3cm)
5	44" – 51" (111.8–129.5cm)
6	50" – 56" (127.0–142.2cm)
EXT.	56" + (142.2cm+)


1


2


3


4


5


6


7


8


9


10

11

12

13


14


15


16


17

18

19

20

21

22

23

OPTIONAL

ENGLISH

INTENDED USE/INDICATIONS:
The intended use of the Ascend and Pinnacle lumbar orthoses is to provide stabilization and motion restriction of the lumbar spine for treatment of trauma related injury, acute post-surgical recovery, and the management of chronic conditions of the spine. Indicated for chronic low back pain, degenerative disc disease, severe DJD, radiculopathy, spondylosis, spinal stenosis, post-operative support, fracture management, post laminectomy syndrome, spondylolisthesis, decompressive procedures, decondensed trunk musculature, sprains/strains, bulging/herniated disc, facet syndrome, IDET procedures, SI dysfunction, sacroiliitis.

INTENDED USERS:
The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use

CONTRAINDICATIONS:
Not Applicable.

⚠ WARNINGS AND PRECAUTIONS:
- Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions.
- If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your licensed healthcare professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.
- This device is offered for sale by or on the order of a licensed health care professional.
- For single patient use only.
- Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
- Follow the licensing healthcare practitioner's protocol regarding method, frequency and duration of administration.

INSTRUCTIONS FOR USE:
MEASUREMENT

1. Measure patient's waist and refer to sizing chart for corresponding size.

ADJUST THE WAIST BELT

Default setting is size 3-4 with belts in the inner slot.

For Sizes 1–3 (Inner Slot):

2. Open the brace and remove the front (anterior) panel.
3. Loosen the black knob and lift it's corner of the mesh cover.
4. Flip the pulley plate over the mesh.
5. Detach the belt from the pulley plate
6. Slide the belt panel through the slot and fold it over until the end aligns with the desired size number. Ensure the belt follows the pulley plate's contour and attach it to the loop.

For Sizes 4–6 (Outer Slot):

7. Remove the hook piece from the pulley plate.
8. Above the belt panel from the inner to the outer slot.
9. Slide and fold the belt to align with the correct size number. Attach the belt end to the loop.
10. Return the pulley plate to its original position and secure the mesh cover.
Repeat the same steps for the other side.
11. Re-attach the belt to the front panel and pull the belt ends away from each other to extend the pulley to its tightest.

LORDOTIC PADS

12. Located inside the back pad. Can be removed or repositioned for extra lower back support

XIPHOID EXTENSION (IF INCLUDED)

13. Align front panel with the lower abdomen.
14. Use the included wrench to adjust the metal plate height.
15. Secure the plate to rest on the xiphoid (lower sternum).

APPLYING THE BRACE

16. Loosen pulley handles are near the mesh cover on the back per #11. (Optional) Attach rigid lateral panels inside the brace.
17. Attach the waist belt to the front panel and secure the hook strap.
18. Align the center of the back pad with the spine and bottom edge over lumbosacral joint.
19. Fasten the brace by pulling the right belt over the left. Adjust belt angle for comfort and fit.
20. Pull handles outward to tighten the brace.
21. Attach pull handles to the belt.

REMOVING THE BRACE

22. Detach pull handles and place them near the back.
23. Undo the front hook-and-loop fastener.

USE AND CARE:
Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or tumble dry.

MATERIAL CONTENTS:
Aimesh®, Aluminum, Cotton, High Density Polyethylene (HDPE), Neoprene, Nylon

WARRANTY:
Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:
Dispose of product according to local regulations.

NOTICE TO USERS: If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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