


Size	Waist Size
1	22" – 28" (55.9–71.1cm)
2	27" – 33" (68.6–83.8cm)
3	32" – 39" (81.3–99.1cm)
4	38" – 45" (96.5–114.3cm)
5	44" – 51" (111.8–129.5cm)
6	50" – 56" (127–142.2cm)
EXT.	56" + (142.2cm+)

ENGLISH

INTENDED USE/INDICATIONS:
The intended use of the Pinnacle Lumbar Orthoses is to provide stabilization and motion restriction for the treatment of human-related injury, acute post-surgical recovery, and the management of chronic conditions of the spine. Indicated for degenerative disc disease, spinal stenosis, spondylosis, bulging/herniated disc, severe DJD, decompressive procedures, kyphosis, scoliosis, kyphosis secondary to osteoporosis, spondylolisthesis, post-operative support and fracture management.

INTENDED USERS:
The intended user should be a licensed professional, the patient, or patient's caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:
Not Applicable.

⚠ WARNINGS AND PRECAUTIONS:
- Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions.
- If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your licensed health care professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.
- Follow application and care process described in instructions.
- For single patient use only.
- Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
- Follow the licensing healthcare practitioner's protocol regarding method, frequency and duration of administration.

INSTRUCTIONS FOR USE:
MEASUREMENT
1. Measure patient's waist and refer to sizing chart for corresponding size.

ADJUSTING THE WAIST BELT SIZE
Default setting is size 3 with belts in the inner slot.
For Sizes 1–3 (Inner Slot):
2. Open brace and remove front panel.
3. Flip it over and lift the mesh cover.
4. Flip the pulley plate.
5. Unhook the belt from the pulley plate.
6. Slide and fold the belt to match the correct size number. **Reattach belt to the loop, following the pulley shape.**
For Sizes 4–6 (Outer Slot):
7. Allow belt from inner to outer slot.
8. Slide and fold belt to match size number. Reattach belt to loop.
9. Put pulley plate back and close mesh cover.
Repeat on the other side.
Lordotic Pads:
10. Lordotic Pads (Lower Back Support)
- Install into the back pad.
- Can be removed or moved for extra support.

ANTERIOR PANEL ADJUSTMENT
11. Hold front panel against abdomen, bottom edge at pubic bone.
12. Pull the pull lever out and slide the metal bar and adjust height.
13. The steroid plate should sit at the sternal notch (top of chest).

APPLYING THE BRACE
Make sure pull handles are near the mesh cover on the back.
14. (Optional) Insert side panels for extra support.
15. Attach waist belt to front panel and secure the strap.
16. Center back pad on spine, bottom edge over lower back.
17. Pull right hand over left and adjust for comfort.
18. Buckle shoulder straps to front panel.
a. Straps can go over shoulders or under arms.
19. Pull handles outward to tighten.
20. Attach handles to the belt.

REMOVING THE BRACE
21. Unhook pull handles and place near the back.
22. a. Unbuckle shoulder straps.
b. Undo front hook-and-loop closure.
You can leave the front panel attached when not wearing.

USE AND CARE:
Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or bleach.

MATERIAL CONTENT:
Airmesh®, Aluminum, Cotton, DIANIT Resin, High Density Polyethylene (HDPE), Neoprene, Nylon, Steel

WARRANTY:
Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:
Dispose of product according to local regulations.

NOTICE TO USERS: If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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