




Size	Waist Size
1	24" – 34" (61.0–86.4cm)
2	33" – 40" (83.8–101.6cm)
3	39" – 46" (99.1–116.8cm)
4	45" – 50" (114.3–127cm)
EXT.	50" + (127cm+)


1


2


3-4


5

6

7

8

9

10

ENGLISH

INTENDED USE/INDICATIONS:

The intended use of these lumbar products is to provide stabilization and motion restriction for the treatment of trauma related injury, acute post-surgical recovery, and the management of chronic conditions of the spine, limited for chronic low back pain, degenerative disc disease, severe DJD, radiculopathy, spondylosis, spinal stenosis, post-operative support, fracture management, post laminectomy syndrome, spondylolisthesis, decompressive procedures, deconditioned trunk musculature, spasm/strains, bulging/herniated disc, facet syndrome, IDET procedures, SI dysfunction, sacroiliitis.

INTENDED USERS:

The intended user should be a licensed professional, the patient or patient's caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:

Not Applicable.

⚠ WARNINGS AND PRECAUTIONS:

- Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions.
- If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reaction, discontinue while using this product, immediately consult your licensed health care professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification can also aid in the overall part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.
- Federal law restricts this device to sale by or on the order of a licensed health care practitioner.
- The frequency and duration of use shall be prescribed by a healthcare professional.
- Follow application and process described in instructions.
- For single patient use only.
- Do not use over open wounds.
- Do not use this device if damaged.

INSTRUCTIONS FOR USE:

MEASUREMENT

1. Measure patient's waist and refer to sizing chart (above) for corresponding size.

ADJUSTING THE WAIST BELT SIZE

Note: The initial configuration is set for size 2.

Open the brace and lay it face down

2. Flip the belt to find sizing numbers

3. Disengage hook and loop fastener and move each end to the desired size.

4. Repeat on both sides

Note: Before trimming any material flip the belt over and wrap around patient's waist to estimate fit, then adjust accordingly.

APPLYING THE BRACE

Ensure pull handles are near the back of the nylon cover

5. Align back pad center with the spine and bottom edge over the lumbosacral joint.

6. Fasten brace by pulling right over left.

7. Adjust tension by pulling handles outward

8. Attach handles to belt.

REMOVING THE BRACE

9. Detach pull handles and place near the back

10. Unfasten front hook and loop to remove.

The brace can be trimmed, bent, molded, assembled, or otherwise customized to fit a specific patient. See product website for further instructions.

USE AND CARE:

Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or tumble dry.

MATERIAL CONTENT:

Cotton, DIANIT resin, High Density Polyethylene (HDPE), Nylon, Polyester, Polyurethane

WARRANTY:

Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:

Dispose of product according to local regulations.

NOTICE TO USERS: If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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