




Size	Waist Size
1	24" – 34" (61.0–86.4cm)
2	33" – 40" (83.8–101.6cm)
3	39" – 46" (99.1–116.8cm)
4	45" – 50" (114.3–127cm)
EXT.	50" + (127cm+)


1


2


3-4


5


6


7


8


10


11


13


15


17


18

12

14

16

9

19

ENGLISH

INTENDED USE/INDICATIONS:
The intended use of these lumbar products is to provide stabilization and motion restriction for the treatment of traumatic and degenerative disorders of the post-surgical recovery, and the management of chronic conditions of the spine. Indicated for chronic low back pain, degenerative disc disease, severe DJD, radiculopathy, spondylosis, spinal stenosis, post-operative support, fracture management, post laminectomy syndrome, spondylolisthesis, decompressive procedures, deconditioned trunk musculature, sprains/strains, bulging/herniated disc, facet syndrome, IDET procedure, kyphosis secondary to osteoporosis, scoliosis with musculoskeletal pain.

INTENDED USERS:
The intended user should be a licensed professional, the patient or patient's caregiver. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

CONTRAINDICATIONS:
Not applicable.

⚠ WARNINGS AND PRECAUTIONS:
- Carefully read fitting instructions and warnings prior to use. To ensure proper performance of the brace, follow all instructions.
- If you experience, or have an increase in pain, swelling, skin irritation, or any adverse reactions while using this product, immediately consult your licensed health care professional.
- This device will not prevent or reduce all injuries. Proper rehabilitation and activity modification are also an essential part of a safe treatment program. Consult your licensed health care professional regarding safe and appropriate activity level while wearing this device.
- Federal law restricts this device to sale by or on the order of a licensed health care practitioner.
- The frequency and duration of use shall be prescribed by a healthcare professional.
- This device is offered for sale by or on the order of a licensed health care professional.
- Follow application and fit process described in instructions.
- For single patient use only.
- Do not use over open wounds.
- Do not use this device if damaged.

INSTRUCTIONS FOR USE:
MEASUREMENT
1. Measure patient's waist and refer to sizing chart for corresponding size.

ADJUSTING THE WAIST BELT SIZE
2. To adjust the waist belt size, open the brace and lay it face down.
3. Flip the waist belt over to find the sizing numbers. Disengage the hook and loop and move one belt end to the desired sizing number.
4. Repeat on the other belt end to complete the sizing process.

CLINICAL TIP
- Flip waist belt over and apply to patient's waist.
- Have patient hold the front of the waist belt to estimate appropriate length for sizing. Filter-adjust belt length for patient's size.
- Flip waist belt back over and apply to patient.

APPLYING THE BRACE
Before applying the brace, make sure the pull handles are positioned on the belt and the sternum strap is unbuckled.

5. Put the brace on by placing the straps over the patient's shoulders.
6. Align the center of the back pad to the patient's spine, and align the bottom edge of the back panel over lumbosacral joint. Refer to picture 6 for proper orientation.
7. Fasten the brace around the patient's waist by pulling the right belt over the left belt. Angle belts up or down to accommodate patient's anatomy.
8. Adjust the tension of the brace by pulling the pull handles away from the body.
9. Attach the buckles to the belt.

10. If no further adjustment is required, buckle the sternum strap and tighten as needed.

ADJUSTING THE POSTERIOR THORACIC EXTENSION AND STRAPS
Make sure the belt is properly fitted before adjusting the posterior thoracic extension (PTE).

11. Unbuckle the sternum strap.
12. While the waist belt is on the patient, place hand between the back pad and PTE to separate the hook attachment.
13. Slide PTE to the desired position and secure in place using hook attachment.
14. Ensure the clips at the base of the shoulder straps are buckled. Adjust the straps as necessary to fit the patient.
15. Re-buckle the sternum strap.
16. To adjust the height of the sternum strap, slide the attachment buckles up or down.

NOTE: Any excess strap length can be secured using the elastic bands.

REMOVING THE BRACE
17. Unbuckle the pull handles from the belt and position them close to the back of the brace.
18. Unbuckle the sternum strap.
19. Detach hook and loop at the front to remove the brace.

Do brace can be trimmed, bent, molded, countersunked, or otherwise customized to fit a specific patient. See product website for further instructions.

USE AND CARE:
Hand wash in 30°C (86°F) water with mild soap, rinse thoroughly. Rinse and lay flat to air dry. Do not machine wash or tumble dry.

MATERIAL CONTENT:
Cotton, High Density Polyethylene (HDPE), Nylon, Polyester, Polyurethane

WARRANTY:
Breg products come with a limited warranty against manufacturing and material defects. No other warranty, expressed or implied, is given and no affirmation of or by seller, by words or action, will constitute a warranty. Please see www.breg.com

DISPOSAL:
Dispose of product according to local regulations.

NOTICE TO USERS: If any serious incident has occurred in relation to the device, report to Breg and the competent authority of the Member State in which the user and/or patient is established.

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