

Lace-Up Ankle with Tibia Strap

Application Instructions



1. Open brace and place foot inside. Tongue should be between foot and laces. Lace and tighten.



2. Place inside strap across top of foot and then under the heel. Attach loop end of strap to hook side of brace.



3. Repeat with the outside strap and attach.



4. Sitting down and with your foot flat on the floor, tighten the straps. To free the straps use the loops at the end of the straps, pull up and then re-attach straps to hook material.



5. When the brace feels comfortable and secure, wrap the elastic top strap over lace and straps. Re-adjust if necessary.

*Note - A white shoe lace has been used to show the application instructions more clearly. However, the Lace-up Ankle with Tibia Strap comes with a black shoe lace.

Measure circumference (with sock on) from heel base around the top of the foot.

Hand Wash • Mild Detergent • Air Dry



Manufactured For:

BREG, Inc. 2611 Commerce Way | Vista, CA 92081

P: 800.321.0607 | F: 800.FAX.BREG

www.breg.com




Lace-Up Ankle with Tibia Strap

Application Instructions



1. Open brace and place foot inside. Tongue should be between foot and laces. Lace and tighten.



2. Place inside strap across top of foot and then under the heel. Attach loop end of strap to hook side of brace.



3. Repeat with the outside strap and attach.



4. Sitting down and with your foot flat on the floor, tighten the straps. To free the straps use the loops at the end of the straps, pull up and then re-attach straps to hook material.



5. When the brace feels comfortable and secure, wrap the elastic top strap over lace and straps. Re-adjust if necessary.

*Note - A white shoe lace has been used to show the application instructions more clearly. However, the Lace-up Ankle with Tibia Strap comes with a black shoe lace.

Measure circumference (with sock on) from heel base around the top of the foot.

Hand Wash • Mild Detergent • Air Dry



Manufactured For:

BREG, Inc. 2611 Commerce Way | Vista, CA 92081

P: 800.321.0607 | F: 800.FAX.BREG

www.breg.com

