

# BACK SUPPORT with POCKET

## FEATURES

- Elastic with 1/16" nylon two sided neoprene pocket to accommodate moldable insert
- Two 4" side pulls overlap at the front for increased abdominal support
- Neutral flesh tone helps maintain low profile
- Allows full-range of motion to help prevent muscle atrophy
- Concentrates heat on lower back
- No stays to irritate while sitting

**INDICATIONS:** Connective tissue and ligament strain where support and heat are prescribed.

**CONTRAINDICATIONS:** If irritations occur, apply talcum powder or under garment between skin and support. If irritation does not clear up, discontinue use. Do not wear while sleeping. Persons with circulatory problems SHOULD NOT wear this product.

**WASHING INSTRUCTIONS:** Join hook and loop closure. Wash in mild detergent or clear, cool water. Air Dry (no dryer). Do not Bleach.

**NEOPRENE WARNING:** Contact of this material with the skin may cause irritation. Persons with known susceptibility to dermatitis or with broken skin should avoid contact with this material. Anyone who develops skin irritation from contact with this material should discontinue its use immediately and consult a physician.



Manufactured For:

Breg, Inc.  
2611 Commerce Way  
Vista, CA 92081 U.S.A.  
P: 800-321-0607  
F: 800-329-2734

www.breg.com  
Breg is an Orthofix company

# BACK SUPPORT with POCKET

## FEATURES

- Elastic with 1/16" nylon two sided neoprene pocket to accommodate moldable insert
- Two 4" side pulls overlap at the front for increased abdominal support
- Neutral flesh tone helps maintain low profile
- Allows full-range of motion to help prevent muscle atrophy
- Concentrates heat on lower back
- No stays to irritate while sitting

**INDICATIONS:** Connective tissue and ligament strain where support and heat are prescribed.

**CONTRAINDICATIONS:** If irritations occur, apply talcum powder or under garment between skin and support. If irritation does not clear up, discontinue use. Do not wear while sleeping. Persons with circulatory problems SHOULD NOT wear this product.

**WASHING INSTRUCTIONS:** Join hook and loop closure. Wash in mild detergent or clear, cool water. Air Dry (no dryer). Do not Bleach.

**NEOPRENE WARNING:** Contact of this material with the skin may cause irritation. Persons with known susceptibility to dermatitis or with broken skin should avoid contact with this material. Anyone who develops skin irritation from contact with this material should discontinue its use immediately and consult a physician.



Manufactured For:

Breg, Inc.  
2611 Commerce Way  
Vista, CA 92081 U.S.A.  
P: 800-321-0607  
F: 800-329-2734

www.breg.com  
Breg is an Orthofix company