
Now
Available in
X-Small
Skyrocket
to Recovery
Designed for little feet, built for
big adventures–helping kids recover
with comfort, style, and durability.
Vectra® Wee Walker
BREG

Introducing Vectra Wee Walker

The Vectra Wee Walker is a pediatric walking boot designed for both comfort and durability. Built to support active children, it combines sturdy engineering to withstand the wear and tear of everyday play. Its adjustable fit ensures a secure, comfortable experience, to help with recovery while staying active. Available with fun, space-themed patches, it's the perfect companion to keep them moving forward in their recovery.

1. **SIDE-MOUNTED AIR PUMP AND RELEASE VALVE** allows effortless adjustments on pneumatic models
2. **LIGHTWEIGHT DESIGN** ensures ease of movement while withstanding the wear and tear of everyday play
3. **HOOK & LOOP STRAPS** for a secure, customizable fit
4. **DURABLE NON-SLIP TREAD** for enhanced traction and stability on various surfaces
5. **ROCKER BOTTOM** promotes natural gait for improved mobility



Vectra Wee Walker Options:

Available in both air and non-air models with size specific designs: the small with a spaceship design, the medium with an astronaut design, and the large with a rocketship design.



SMALL AIR AND NON-AIR



MEDIUM AIR AND NON-AIR




LARGE AIR AND NON-AIR

Vectra Wee Walker Patches:

Add a little extra fun with optional patches! Choose from the pink pack of 3 or the gray pack of 3, each featuring fun space-themed designs to give each walker boot a personalized touch.



PINK PATCHES



GRAY PATCHES



Breg, Inc, 2382 Faraday Ave, Ste # 300,
Carlsbad, CA 92008, USA
P: 800.321.0607, F: 760.795.5295,
E: orderprocess@BREG.com, W: BREG.com

VECTRA WEE WALKER BOOTS

BL0100-010	**X-Small** Vectra Wee Walker Boot
BL0100-020	**Small** Vectra Wee Walker Boot
BL0100-030	**Medium** Vectra Wee Walker Boot
BL0100-040	**Large** Vectra Wee Walker Boot
BL0200-020	**AIR Small** Vectra Wee Walker Boot
BL0200-030	**AIR Medium** Vectra Wee Walker Boot
BL0200-040	**AIR Large** Vectra Wee Walker Boot

REPLACEMENT KITS

BL0006-010	**X-Small** Replacement Kit
BL0006-020	**Small** Replacement Kit
BL0006-030	**Medium** Replacement Kit
BL0006-040	**Large** Replacement Kit
BL0007-020	**AIR Small** Replacement Kit
BL0007-030	**AIR Medium** Replacement Kit
BL0007-040	**AIR Large** Replacement Kit

PATCH ACCESSORIES

BL0000-100	Pink Patches - Set of 3
BL0000-200	Gray Patches - Set of 3

SIZE CHART

BOOT SIZE	FOOT LENGTH
XS	Up to 5 inches (12.7 cm)
S	5 - 6 inches (12.7 - 15.2 cm)
M	6 - 7 inches (15.2 - 17.8 cm)
L	7 - 8 inches (17.8 - 20.3 cm)

©2026 BREG, Inc. All rights reserved. Vectra is a registered trademark of BREG, Inc.
P/N: 10.000019 Rev C 07/26

CE MD EC REP

MDSS GmbH
Schiffigraben 41
D-30175 Hannover
Germany